



FACE SHIELD

mVISOR7001

technical specification



TECHNICAL SPECIFICATION: mVISOR 7001

Use

mVisor 7001 is a protective face shield that offers approved protection against droplets and liquid splashes. It should be used in environments where the user is exposed to liquid splashes and/or droplets.

The product should only be used in accordance with the user manual.

- All-over protection for your face against splashes of liquid and droplets
- Good visibility as the face shield material is clear, transparent and prevents fogging
- Padded with foam and with adjustable straps for an optimal fit and high comfort
- Low Weight: <200 g

Product description


Name of product: mVisor 7001
Product description:
Protective face shield

Labelling

Face shield: Marka Pac 1 S
Head strap: Marka Pac EN 166 3 S

EC declaration of conformity
can be downloaded at
mshield.se/docs

According to requirements

Regulation (EU) 2016/425	•
CE EN 166:2001 	•
No. 20/3121/00/0161 (AITEX)	•
Plaza Emilio Sala No. 1	•
Alcoi, Spain, Notified Body 0161	•

Teknisk specifikation

Lens	Adjustable head strap	Foam
Thickness: 1,8 mm	Thickness: 1,2 mm	Thickness: 30 mm
Length: adjustable 219 mm	Width: 25 mm	Width: 30 mm
Weight: 146 gr	Weight: 24,2 gr	Length: 200–330 mm
Material: polycarbonate, PC	Material: polypropylene, PP	Weight: 6,2 gr
		Material: BPI PolyEther+PolyEster

Care instructions

mVisor 7001 can be used for as long as the shield provides good visibility and is undamaged. Regularly check that the head strap and lens are free from cracks.

- Keep in sealed original packaging for best durability
- Do not leave the face shield in direct sunlight
- Damaged or scratched shields must be replaced

Cleaning instructions

1. Disinfect your hands with hand sanitizer
2. Put on protective gloves
3. Disinfect and clean the shield and head strap with soap and water or use a surface disinfectant
4. Hang up or put the face shield in a place where it can dry, or wipe using a surgically clean disposable wipe
5. When dried and disinfected, put the face shield in a clean plastic bag or original packaging before reuse



Recycling

At the end of the face shield's useful life, disassemble in reverse order to the assembly instructions. Shield and strap can be sorted at source as rigid plastic waste in accordance with local recycling rules.



Documentation available for download

Please visit: mshield.se/docs

Contact

mShield | Anders Anderssons väg 12, 285 35 | Markaryd, Sverige | info@mshield.se | +46 433 128 45